

诺禾致源数据结果说明文档 ---原始数据结果



高通量测序(如 Illumina HiSeq XTen/2500/2000/MiSeq 等测序平台)得到的原始图像数据文件经碱基识别(Base Calling)分析转化为原始测序序列(Sequenced Reads),我们称之为Raw Data 或 Raw Reads

本文件夹下包含每个样本(文库名)的原始数据结果文件:

1.后缀 fq.gz 文件 高通量测序的原始测序序列压缩文件

2.MD5.txt fq.gz 的 MD5 值,用于检查文件的完整性

文件名中 1,2 分别代表 Paired-end 两端序列文件

结果文件说明

1.fq.gz 为高通量测序的原始测序序列,结果以 FASTQ 文件格式存储。包含测序序列的序列信息和对应的测序质量信息。FASTQ 文件中每个 read 由四行描述。其格式如下:

@HWI-ST1276:71:C1162ACXX:1:1101:1208:2458 2:N:0:CGATGT CTGGCTCCGGAGGGGATGGAGGCGGCACTCCCGCCAAGGATGCGTTGGGAAACGACG TCGTTGCAGTCGAATGGCTCAAAACACACGGGCCCGGTGACCGG

+

BCBFFFFDHHHHHJJ?EAGIIIAHIJIIGHHHBEDCDDD;>>BD?BDAD<><?BDB@5<BBDCDDDCDDBDCCDDD8?AAB9>B55>BB5904@BB

其中第一行以"@"开头, 随后为 Illumina 测序标识符 (Sequence Identifiers) 和描述文字 (选择性部分);

第二行是碱基序列;

第三行以"+"开头,随后为 Illumina 测序标识别符(选择性部分);

第四行是对应碱基的测序质量,该行中每个字符对应的 ASCII 值减去 33,即为对应第二行碱基的测序质量值。

Illumina 测序标识别符(Sequence Identifiers)详细信息如下:

HWI-ST1276	Unique instrument name
71	D ID

71 Run ID
C1162ACXX Flowcell ID
1 Flowcell lane

Tile number within the flowcell lane
1208 'x'-coordinate of the cluster within the tile
2458 'y'-coordinate of the cluster within the tile

2 Member of a pair, 1 or 2 (paired-end or mate-pair reads only)

N Y if the read fails filter (read is bad), N otherwise

0 when none of the control bits are on, otherwise it is an even number

CGATGT Index sequence